

MAC100

ARCATRON ASSISTED COMMODE WHEELCHAIR [MAC100]



MAC100 is an attendant propelled commode wheelchair designed for caregivers to enhance hygiene routine for individuals with limited mobility. It allows access to toilets in a secured and convenient manner. Comfortable padded toilet seats and front open cut out makes cleaning access very easy and practical. Automotive grade coated steel frame for increased corrosion resistance. A very special feature allows it to be upgraded to self propelled wheelchair by adding Arcatron MSC100 expansion pack.

SALIENT FEATURES



Can Roll Over Commodes Upto 17" High



Automotive Grade Powder Coating: Water Resistant Frame, Stainless Hardware.



Swing Back Swivel Armrests for Easy Bed Transfer



Removable Toilet Pan (*optional accessory)



Swing Back Footrests for Ease of Transfer and Easy Storage



Toolfree Removable Backrest for Easy Storage



Medical Grade Corrosion Resistant Nylon Casters with Brakes and Swivel Lock. Also prevents Hair Entry

BEST IN-CLASS UPHOLSTERY

Comfortable, Soft-Padded, Anti Bacterial and Waterproof - Silicon Coated Closed Cell Molded PU Cushions



Keyhole Seat with Front Cutout for Easy Cleaning Access



Front Cutout Cover for Hassle Free Transfer (*optional accessory)



Keyhole Seat Cover to Use the Chair as Regular Wheelchair (*optional accessory)



Tension Adjustable High Neck Support Nylon Backrest

ARCATRON ASSISTED COMMODE WHEELCHAIR(MAC100)

WHAT DO YOU GET??

BEST IN-CLASS FEATURES.

Compare the Arcatron MAC100 with other available commode chairs in the same category.



Products	ARCATRON MAC100	COMPETITOR MODEL A	COMPETITOR MODEL B	COMPETITOR MODEL C
Features				
Comfortable and padded toilet seat	✓	✗	✗	✗
Front open seat for easy cleaning access	✓	✗	✗	✗
Tension adjustable high neck backrest for comfort	✓	✗	✗	✗
Swing away armrests for easy transfers	✓	✗	✗	✓
Medical Grade caster wheels for corrosion resistance/hair entry	✓	✗	✗	✗
Upgradeable to Self Propelled Wheelchair ***	✓	✗	✗	✗

*** Arcatron MAC100 assisted wheelchair can be upgraded to Arcatron MSC100 Self Propelled Wheelchair without any hassle. Just add the MSC100 Expansion Pack. Live a more independent life.



SPECIFICATIONS



Seat Height	21"
Seat Width	19"
Seat Depth	16.5"
Seat Recline Posterior Angle	9 degrees
Overall Width	23"
Overall Height	42"
Overall Length	28"
Backrest Width	19"
Backrest Height	22"
Compatible Commode Height	16"
Compatible Commode Width (Max)	16"
Armrest Legth	19"
Weight Capacity	90 kgs
Front Caster Wheels (2No.s)	5" Dia (without brake and swivel lock)
Rear Caster Wheels (2No.s)	5" Dia (with brake and swivel lock)

**specifications are subject to change without notice

contact no.: +91-8600054430

website: www.arcatron.com