


# SSS100

## SPECIFICATIONS (Self-Propelled)

PARAMETERS	NOTATION	MODEL	
		SSS100	SSS100N (Without Magwheel Pushrims)
Dimension	H* x W x D	34-38" x 26" x 42"	34-38" x 24-25.5" x 42"
Weight	-	22.75 Kg	22.75 Kg
Weight Capacity	-	120 Kg	120 Kg
Seat Dimensions	C x A	19" x 16"	19" x 16"
Width Between Arms	E	16.5"	16.5"
Commode Clearance	M* x G	(17"-21") x 16.5"	(17"-21") x 16.5"
Seat Depth	A	16"	16"
Foot Plate Distance	J	19-21"	19-21"
Seat to Arm Height	F	8"	8"
Seat Height	B*	20" - 24"	20" - 24"
Back Dimensions	K x L	17" x 16"	17" x 16"



\*For 19" Commode clearance



# SAS100

## SPECIFICATIONS (Attendant-Propelled)

PARAMETERS	NOTATION	MODEL	
		SAS100	SAS100 (With Magwheel Breaks)
Dimension	H* x W x D	34-40" x 20.5" x 39"	34-40" x <u>24"</u> x 39"
Weight	-	19.75 Kg	19.75 Kg
Weight Capacity	-	120 Kg	120 Kg
Seat Dimensions	C x A	19" x 16"	19" x 16"
Width Between Arms	E	16.5"	16.5"
Commode Clearance	M* x G	(16"-22") x 16.5"	(16"-22") x 16.5"
Seat Depth	A	16"	16"
Foot Plate Distance	J	19"-21"	19"-21"
Seat to Arm Height	F	8"	8"
Seat Height	B*	19"-25"	19"-25"
Back Dimensions	K x L	17" x 16"	17" x 16"

